KAGOSHIMA CASE REPORT 02

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1. Kaneko M, Katanuma A, Maguchi H et al. Prospective, randomized, comparative study of delineation capability of radial scanning and curved linear array endoscopic ultrasound for the pancreaticobiliary region. *Endosc Int Open* 2014; 02: E160-170.
2. Tanaka M, Castillo CF, Kamisawa T et al. Revisions of international consensus Fukuoka guidelines for the management of IPMN of the Pancreas. *Pancreatology* 2017; 17: 738-753.
3. Kamata K, Kitano M, Omoto S et al. Contrast-enhanced harmonic endoscopic ultrasonography for differential diagnosis of pancreatic cysts. *Endoscopy* 2016; 48: 35-41.
4. Tanno S, Nakano Y, Koizumi K et al. Pancreatic ductal adenocarcinomas in long-term follow-up patients with branch duct intraductal papillary mucinous neoplasms. *Pancreas* 2010; 39: 36-40.
5. Kamata K, Kitano M, Kudo M et al. Value of EUS in early detection of pancreatic ductal adenocarcinomas in patients with intraductal papillary mucinous neoplasms. *Endoscopy* 2014; 46: 22-29.
6. Tanaka M, Castillo CF, Adsay V et al. International consensus guidelines 2012 for the management of IPMN and MCN of the Pancreas. *Pancreatology* 2012; 12: 183-197.
7. Kimura W, Moriya T, Hirai I et al. Multicenter study of serous cystic neoplasm of the Japan pancreas society. *Pancreas* 2012; 41: 380-387.
8. 花田敬士, 南智之, 清水晃典, ほか: 精査による臨床像の差違—Solid Pseudopapillary Neoplasm (SPN)—. *胆と膵* 2018; 39: 547-51.

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1. Vilmann P, Jacobsen GK, Henriksen FW, et al. Endoscopic ultrasonography with guided fine needle aspiration biopsy in pancreatic disease. *Gastrointest Endosc.* 1992; 38: 172-173.
2. 抗血栓薬服用者に対する消化器内視鏡診療ガイドライン 直接経口抗凝固薬 (DOAC) を含めた抗凝固薬に関する追補2017. *日本消化器病学会雑誌* 2017; 59: 1547-1558.
3. Bang JY, Hawes R, Varadarajulu S. A meta-analysis comparing ProCore and standard fine-needle aspiration needles for endoscopic ultrasound-guided tissue acquisition. *Endoscopy* 2016; 48: 339-349.
4. Facciorusso A, Wani S, Triantafyllou K, et al. Comparative accuracy of needle sizes and designs for EUS tissue sampling of solid pancreatic masses: a network meta-analysis. *Gastrointest Endosc* 2019; 90: 893-903.e897
5. Han S, Bhullar F, Alaber O, et al. Comparative diagnostic accuracy of EUS needles in solid pancreatic masses: a network meta-analysis. *Endosc Int Open.* 2021; 9: E853-E862.
6. Gimeno-Garcia AZ, Paquin SC, Garipey G et al. Comparison of endoscopic ultrasonography-guided fine-needle aspiration cytology results with and without the stylet in 3364 cases. *Dig Endosc* 2013; 25: 303-307.
7. Nakai Y, Isayama H, Chang KJ, et al. Slow pull versus suction in endoscopic ultrasound-guided fine-needle aspiration of pancreatic solid masses. *Dig Dis Sci.* 2014; 59: 1578-85.
8. Larghi A, Noffsinger A, Dye CE et al. EUS - guided fine needle tissue acquisition by using high negative pressure suction for the evaluation of solid masses: a pilot study. *Gastrointest Endosc.* 2005; 62: 768-774.
9. 山雄健次. 超音波内視鏡ガイド下穿刺術. 私のコツ. *消化器画像* 2007; 9: 98-104.
10. 入澤篤志. 超音波内視鏡ガイド下穿刺術. 私のコツ. *消化器画像* 2007; 9: 199-205.
11. Bang JY, Magee SH, Ramesh J et al. Randomized trial comparing fanning with standard technique for endoscopic ultrasound-guided fine-needle aspiration of solid pancreatic mass lesions. *Endoscopy* 2013; 45: 445-50.
12. Iglesias-Garcia J, Dominguez - Munoz JE, Abdulkader I, et al. Influence of on-site cytopathology evaluation on the diagnostic accuracy of endoscopic ultrasound-guided fine needle aspiration (EUS-FNA) of solid pancreatic masses. *Am J Gastroenterol* 2011; 106: 1705-1710.
13. Crinò SF, Di Mitri R, Nguyen NQ, et al. Endoscopic Ultrasound-guided Fine-needle Biopsy With or Without Rapid On-site Evaluation for Diagnosis of Solid Pancreatic Lesions: A Randomized Controlled Non-Inferiority Trial. *Gastroenterology.* 2021; 161: 899-909
14. Suzuki R, Irisawa A, Bhutani MS et al. Prospective evaluation of the optimal number of 25-gauge needle passes for endoscopic ultrasound-guided fine-needle aspiration biopsy of solid pancreatic lesions in the absence of an onsite cytopathologist. *Dig Endosc.* 2012; 24: 452-456
15. Wang KX, Ben QW, Jin ZD, et al. Assessment of morbidity and mortality associated with EUS-guided FNA: a systematic review. *Gastrointest Endosc.* 2011; 73: 283-290.
16. Owens CL, Peterson D, et al. Effects of transitioning from conventional methods to liquid-based methods on unsatisfactory Papanicolaou tests: results from a multicenter US study. *Cancer Cytopathol.* 2013; 121: 568-75.
17. Hashimoto S, Taguchi H, Higashi M, et al. Diagnostic efficacy of liquid-based cytology for solid pancreatic lesion samples obtained with endoscopic ultrasound-guided fine-needle aspiration: Propensity score-matched analysis. *Dig Endosc.* 2017; 29: 608-616.
18. 日本臨床腫瘍学会・日本癌治療学会・日本癌学会 3学会合同ゲノム医療推進タスクフォース「血中循環腫瘍DNAを用いたがんゲノムプロファイリング検査の適正使用に関する政策提言」(令和3年1月20日) <https://www.jsmo.or.jp/file/dl/newsj/2765.pdf>
19. Bang JY, Hebert-Magee S, Navaneethan U, et al. Randomized trial comparing the Franseen and Fork-tip needles for EUS-guided fine-needle biopsy sampling of solid pancreatic mass lesions. *Gastrointest Endosc.* 2018; 87: 1432-1438.

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1. M. Giovannini, V. Moutardier, C. Pesenti, E. Bories, B. Lelong, and J. R. Delperio, "Endoscopic ultrasound-guided bilioduodenal anastomosis: A new technique for biliary drainage," *Endoscopy*, vol. 33, no. 10, 2001.
2. E. Burmester, J. Niehaus, T. Leineweber, and T. Huetteroth, "EUS-cholangio-drainage of the bile duct: Report of 4 cases," *Gastrointest. Endosc.*, vol. 57, no. 2, 2003.
3. B. A. Holt et al., "Biliary drainage: Role of EUS guidance," *Gastrointest. Endosc.*, vol. 83, no. 1, pp. 160–165, 2016.
4. 伊佐山 浩通, 中井 陽介, 糸井 隆夫, 安田 一朗, 河上 洋, 良沢 昭銘, 北野 雅之, 入澤 篤志, 瀧沼 朗生, 原 和生, 岩下 拓司, 藤田 直孝, 山雄 健次, 吉田 雅博, 乾 和郎, "超音波内視鏡下胆道ドレナージの安全施行への診療ガイドライン：2018," *胆道*, vol. 33, no. 5, pp. 793–816, 2019.
5. S. Mallery, J. Matlock, and M. L. Freeman, "EUS-guided rendezvous drainage of obstructed biliary and pancreatic ducts: Report of 6 cases," *Gastrointest. Endosc.*, vol. 59, no. 1, 2004.
6. N. Fujita et al., "Temporary endosonography-guided biliary drainage for transgastrointestinal deployment of a self-expandable metallic stent," *J. Gastroenterol.*, vol. 43, no. 8, 2008.
7. S. J. Park et al., "Expanding indication: EUS-guided hepaticoduodenostomy for isolated right intrahepatic duct obstruction (with video)," *Gastrointest. Endosc.*, vol. 78, no. 2, 2013.
8. 石渡裕俊, "超音波内視鏡下胆道ドレナージの進歩と課題," *胆道*, vol. 34, no. 1, pp. 76–88, 2020.
9. Hashimoto S, et al. Comparison of recurrent biliary obstruction with the use of metal and plastic stents in endoscopic ultrasonography-guided biliary drainage: A propensity score-matched analysis. *Endosc Ultrasound*. In press

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