

ES DILATOR -Super Soft type-

XIEM
EX

Universal catheter dilator



ES Dilator
Effective Dilator for Interventional EUS

ZEON MEDICAL INC.

Universal catheter dilator



Features

- By making the tubing material more flexible into fluoroethylene, it is softer and smoother to insert than ES dilators (soft types).
- An Imaging ring located 10 mm from the tip can be used to identify the tip position under fluoroscopy. The imaging ring is inside the catheter, so the outside of the catheter is smooth.
- The identification code is printed on the hub (Super soft type only).



SPEC

P / N	Catheter outer diameter	Catheter effective length	Applicable guide wire
DC7F180M(E)	2.45mm(7F)	1800mm	0.64mm(0.025 inch)

attention

- Before using product, be sure to read the attached IFU.
- Product is a custom-made product. For details, please contact your sales representative.

manufacturer

ZEON MEDICAL INC.

[H O] Shin-Marunouchi Center Building, 1-6-2 Marunouchi, Chiyoda-ku, Tokyo 100-0005, Japan
Phone: +81-3-3216-1265
URL : <https://www.zeonmedical.co.jp/e/>

XEMEX is a house mark of Zeon Medical Inc. and is registered in Japan and other countries.