

S-O clip


Traction Clip

XIEM
EX



s-o clip

ZEON MEDICAL INC.




ESD(Endoscopic Submucosal Dissection) Towards securing the visual field and reducing the procedure time.

FEATURE

- It is possible to grasp the mucous membrane with a simple operation that follows the Zeoclip™ technology.
- Continuous traction of the lesion by the spring effect makes it easier to secure the field of view on the dissection surface.
- Because the lesion is constantly traction, the latent maneuver can be reduced, and ESD can be performed safely and in a short time.

HOW TO USE



SPEC

P/N	Spring length	Loop diameter	Packaging unit
TC1H05 (E)	5mm	4mm	10pcs


attention

- Product is used by physicians or physician-supervised healthcare professionals and should be used with adequate training on endoscopic procedures.
- Before using product, be sure to read the attached IFU.

Combined device

ZEACLIP™

ZEACLIP
QUICK, STABLE, EASY
Effective Solution



manufacturer

ZEON MEDICAL INC.

[HO] Shin-Marunouchi Center Building, 1-6-2 Marunouchi, Chiyoda-ku, Tokyo 100-0005, Japan
Phone: +81-3-3216-1265
URL: zeonmedical.co.jp/e/

XEMEX is a house mark of Zeon Medical Inc. and is registered in Japan and other countries.